

Post tagmentation

5 clean up



ST2 TWB	25	Vortex	5 OT-REAGENT OT-REAGENT REAGENT PLATE TRANSFER to OT-PREP	Add water to TM Transfer TM to Reagent Plate	P300 single P300 single	20 ul per sample 90 ul per well
	4	Vortex		6 OT-PREP1 7 OT-PREP1 Transfer from Reagent Plate to TAG1 and mix Thermal cycle	P300 multi Thermal Cycler	30 ul per sample
TAG1 TRANSFER FROM THERMAL CYCLER TO MAGNETIC MODULE						

Reservoir for TWB

1 OT-REAGENT REAGENT PLATE TRANSFER to OT-PREP	Distribute ST2 to Reagent plate	P300 single	10 ul per sample
	2 OT-PREP1 3 OT-PREP1 4 OT-PREP1 5 OT-PREP1 6 OT-PREP1 7 OT-PREP1 Trasfer from Reagent plate to TAG1 and mix Incubate Activate magnetic module Wait for liquid to be clear Discard the supernatant Disable magnetic stand	P20 multi P300 multi	10 ul per sample 5 min 3 min 50 ul
REAGENT PLATE TRANSFER to OT-PREP			
12-WELL RESERVOIR TRANSFER to OT-PREP			
8 9 10 11 12 13	Distribute TWB to 12-well Reservoir Add TWB to TAG1 Mix Engage magnetic module Discard the supernatant Add TWB to TAG1 Mix	P1000/P300 single P300 multi P300 multi P300 multi P300 multi P300 multi	200 ul per sample 100 ul per sample 100 ul 100 ul 100 ul per sample 100 ul

Amplify Tagmented

6 Amplicons

6 Amplicons	EPM	-20	ice	New 15ml tube (MM)
	Index adapters water	-20	Thaw ambient; vortex;	
				
				

Pool and clean up

7 libraries

7 libraries	ITB	25	Vortex before each use	New 1.7ml tube (Pooled ITB) New 200ul tube (Eluted) EtOH Reservoir
	RSB	4	Wait 30 min to bring at	
	EtOH	0	Prepare fresh; 80%	

POOL STRIP TRANSFER TO OT-PREP2 (Magnetic module)					15 ul per well
4 OT-PREP2		Transfer from Pool Strip to Pooled ITB and mix		P300 single	
QUANTIFICATION PLATE TRANSFER TO OT-PREP2					
5 OT-PREP2	Mix ITB			P300 single	
6 OT-PREP2	Add ITB to Pooled ITB and mix			108 ul	
7 OT-PREP2	Incubate			5 min	
9 OT-PREP2	Activate Magnetic Module				
OT-PREP2	Wait for liquid to be clear			5 min	
10 OT-PREP2	Discard supernatant			400 ul	
12 OT-PREP2	Wash beads: add EtOH 80% to Pooled ITB			1000 ul	
13 OT-PREP2	Wait			30 s	
14 OT-PREP2	Discard supernatant			1000 ul	
15 OT-PREP2	Wash beads: add EtOH 80% to Pooled ITB			1000 ul	
16 OT-PREP2	Wait			30 s	
17 OT-PREP2	Discard supernatant			1000 ul	
18 OT-PREP2	Discard supernatant with 20ul pip			P300 single	
OT-PREP2	Disengage magnetic module			P20 single	
19 OT-PREP2	Add RSB to Pooled ITB and mix			P300 single	
20 OT-PREP2	Incubate			55 ul	
21 OT-PREP2	Engage magnetic module			2 min	
22	Wait for liquid to be clear			2 min	
23	Transfer to Eluted tube			50 ul	
P300 single					

Safe stop

Quantify and					
8	Normalize Library	dsDNA HS Reagent	4	Bring at ambient	
		dsDNA HS Buffer	<=30	Bring at ambient	
	See Qubit dsDNA HS	dsDNA HS Std1	4	Bring at ambient	
	Assay Kit User Guide	dsDNA HS Std2	4	Bring at ambient	
				Wait 30 min to bring at room temperature;	
		RSB	4	vortex and invert	
New 3 Qubit Assay Tubes (STD1, STD2, SAMPLE)					
New 1.7ml tube (working solution)					
[New microcentrifuge tube (Diluted)]					
1	OT-REAGENT	Add dsDNA HS Buffer to working solution		P300 single	600 ul circa
		Add dsDNA HS Reagent to working solution and mix			
2	OT-REAGENT			P20 single	1:200
3	OT-REAGENT	Transfer working solution to STD1, STD2		P300 single	rispetto a HS Bu
4	OT-REAGENT	Transfer working solution to SAMPLE		P300 single	190 ul
5	OT-REAGENT	Add dsDNA HS Std1 to STD1		P20 single	198 ul
6	OT-REAGENT	Add dsDNA HS Std2 to STD2		P20 single	10 ul
				P20 single	10 ul
QUANTIFICATION PLATE TRANSFER TO OT-PREP2					
7	OT-PREP2	Add Eluted to SAMPLE		P20 single	2 ul
8	OT-PREP2	Mix STD1		P300 single	
9	OT-PREP2	Mix STD2		P300 single	
10	OT-PREP2	Mix SAMPLE		P300 single	
USER INTERVENTION: Qubit quantification (using STD1, STD2, SAMPLE)					
Result: SAMPLE concentration in ng/ul					
11		Calculate volume for dilution to 4nM			
12	OT-PREP2	Add RSB to Eluted as per dilution		P300 single/P20 single	? ul

Dilute to final loading

9 concentration

See Miniseq Denature and Dilute Libraries Guide

Thaw ambient; then keep 4°C until use; vortex before use
Prepare fresh; 1ml 0.1N

-20

HT1

NaOH water
200mM Tris-HCl pH 7.0

(Questions) RSB ?

Thaw ambient; then keep 4°C until use

-20

New 3 1.7ml tube (NaOH, LIBRARY, LIBDILUTED)
New microcentrifuge tube (TUBE1)

-> Prepare NaOH 0.1N

1 OT-REAGENT
2 OT-REAGENT

P300 single
P300 single

900 ul
100 ul

-> Dilute library to 1nM

3 OT-PREP2
4 OT-PREP2
5
6
7

P20 single
P300 single
P20 single
P20 single

25 ul
75 ul
5 ul
5 ul
5 min

-> Denature library

DENATURE PLATE TRANSFER TO OT-PREP2

Transfer Library to TUBE1
Transfer RSB to TUBE1 and mix
Transfer from TUBE1 to LIBRARY
Transfer from NaOH to LIBRARY and mix
Incubate at room temperature

-> Dilute library to final loading concentration of 1.2pM

8
9
10

P20 single
P300 single
P300 single

5 ul
985 ul
380 ul

-> MiniSeq loading volume is 500ul

11

P300 single

120 ul

Preparation finished -> Load to MiniSeq